

RESOLUTIONS. Intra Spectral Imaging System.

Our ref: ISIS0099 Date: 18/11/2004 Time: 17:27:52



ISIS Image from www.resolutions.org.uk

Test page to show ISIS Image Print.

Any text here will be displayed along side of the image in the report.

You can CUT / COPY and PASTE by using your right mouse button. Word wrap is enabled by default, so only use carriage return when you want to start a new paragraph.

You are limited to approximately 900 characters inc. spaces.



ISIS Image from www.resolutions.org.uk

Test page to show ISIS Image Print.

Any text here will be displayed along side of the image in the report.

You can CUT / COPY and PASTE by using your right mouse button. Word wrap is enabled by default, so only use carriage return when you want to start a new paragraph.

You are limited to approximately 900 characters inc. spaces.

Therapist ref: 1234

Client Name: Niki

This Scan No: 0200

Original Scan No: 0100

Image1 Date: 16/11/2004

Filter used: 003

Image2 Date: 16/11/2004

Filter used: 001

Notes: Test page to show ISIS Image Print.

Any text here will be displayed at the bottom of the report. You can CUT / COPY and PASTE by using your right mouse button. Word wrap is enabled by default, so only use carriage return when you want to start a new paragraph. You are limited to approximately 700 characters inc. spaces.